





Digital Slide Scanners



Application Module

Vision Slide



Object analysis

Objects classification by specified criteria.

Count and classification

Automatic or manual object identification on the digital slide. Measurements of linear and optical parameters for particles.

Marks and comments

Make your comments and put marks on a digital slide.

Administrative Module

Vision Remote



Remote workstation.

Consulting and Education Module

Vision Suite

The screenshot shows the Vision Suite web interface. The browser address bar displays "vision-suite.com". The page shows a list of histology images with the following columns: Image, List, Sample Date, ID, Sample, Size, and Note. The data rows are as follows:

Image	List	Sample Date	ID	Sample	Size	Note
		12/05/2016, 9:49 PM	2850	Histo	71 KB	
		12/05/2016, 9:49 PM	2849	Histo	437.4 KB	
		12/05/2016, 9:44 PM	2848	Histo	345.3 KB	
		12/05/2016, 9:39 PM	2847	Histo	400.3 KB	
		12/05/2016, 7:33 PM	4814	Histo	1.8 KB	
		12/05/2016, 7:33 PM	4858	Histo	905.7 KB	
		12/05/2016, 7:33 PM	4857	Histo	1.8 KB	

Cloud/server for telemedicine and remote consultations with colleagues.

Specifications



Vision Assist
Cell Imaging Analyzer



Vision Pro
Cell Imaging Analyzer



Vision Ultimate
Cell Imaging Analyzer

Application module: Vision Slide	Application module: Vision Slide	Application module: Vision Slide
Working modes: queue (only 4 slides version)	Working modes: queue, random access	Working modes: queue, continuous and random access, STAT testing, 24/7
Automatic scanning	Automatic scanning	Automatic scanning
1 or 4 slides	4 or 8 slides	Up to 200 slides
Manual slide handling	2 cassettes with slides	Automatic slide handling
Microscope for scanning	Microscope for scanning	Microscope for scanning
Personal computer	Personal computer	Personal computer
Monitor	Monitor	Monitor
—	—	Touch-screen monitor for control
—	Built-in barcode reader (optional)	Built-in barcode reader
—	Automatic oil dispenser (optional)	Automatic oil dispenser (optional)
Optical system: 20x*	Optical system: 20x*	Optical system: 10x, 20x*
Bright field	Bright field	Bright field
Köhler illumination, LED	Köhler illumination, LED	Köhler illumination, LED
Bidirectional LIS, LIS2-A2 (ASTM), HL7, Ethernet	Bidirectional LIS, LIS2-A2 (ASTM), HL7, Ethernet	Bidirectional LIS, LIS2-A2 (ASTM), HL7, Ethernet
Art. N.: 71150.00.011 (1 slide) Art. N.: 71450.00.011 (4 slides)	Art. N.: 72452.00.011 (4 slides) Art. N.: 72852.00.011 (8 slides)	Art. N.: 73051.00.011 (200 slides)

* Ask for a full list of available options when ordering

We reserve the right to change specification without notice.



Manufacturer:
West Medica Produktions- und Handels-GmbH
Brown-Boveri-Straße 6, B17-1, 2351 Wiener Neudorf, Austria
tel.: +43 (0) 2236 892465, fax: +43 (0) 2236 892464
vienna@westmedica.com, www.westmedica.com

For distributor